Movora – North American Headquarters USA  
310 Commerce Lake Dr, Ste. 107  
St. Augustine, Florida 32095 USA  
www.movora.com

Movora – European Headquarters Switzerland  
KYON AG  
Hardturmstrasse 103  
8005 Zürich Switzerland  
www.movora.com

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## ENGLISH

### 1. Device Description

The OrthoSta System consists of two sizes of woven Ultra-High Molecular Weight Polyethylene (UHMWPE) tape coated in Healix, a silver ion antimicrobial, prefixed with a stainless-steel blunt needle and titanium suture button. The OrthoSta system is used in conjunction with titanium alloy interference screws, stainless-steel taps, and surgical accessories including a screwdriver and drill bits.

### 2. Indications

The OrthoSta System is intended for veterinary use in the management of ligament and tendon injuries. Applications include stabilization of cranial cruciate injuries, ligament injuries, collateral ligament injuries, common calcaneal injuries, patellar tendon injuries, triceps tendon injuries, hip luxation, and shoulder instability.

### 3. Contraindications

Contraindications include insufficient quality or quantity of bone, blood supply issues or previous infections that will prevent healing, foreign body sensitivity, infection, and or conditions that may prevent the patient from following post op instructions. Do not use for surgeries other than those indicated.

### 4. Adverse Effects

Infection, foreign body reactions, soft tissue inflammation, aseptic loosening, and failure/breakage of the implants may occur.

### 5. Warnings

The OrthoSta System is intended to be used by a trained medical professional and is for veterinary use only. This is a single use device and must not be re-used or re-sterilized. Failure to use the proper drill bit diameters or taps may cause failure or breakage of the implant. The surgeon must be aware that if the implant needs to be removed it will require a second surgery. Explanted devices should be safely disposed of and not re-used. Knowledge of surgical techniques and proper selection and placement of the implants are important considerations in the successful use of the OrthoSta System.

### 6. MRI Safety Information

The OrthoSta interference screws and suture buttons have not been evaluated for safety and compatibility in the MR environment. Suture manufactured from UHMWPE is MR safe.

### 7. Precautions

Veterinarian medical professionals are advised to review the technical brochure prior to surgery. Surgeons may apply their professional judgment in the selection of the appropriate implants. Ensure the indicated drill bits and taps are used for the corresponding size tape and interference screw. Respect the weight limits advised in the technical brochure.

### 8. Package and Labeling

Devices should only be accepted if the packaging and labeling arrive intact. Please contact customer service if the package has been opened or altered upon arrival. Do not use any device if the package or implant is damaged.

### 9. Sterilization

The OrthoSta tape is provided sterile using Ethylene Oxide. Check the labeling for more information and expiration date. This device is single use only; do not re-sterilize. Certain instruments and accessories may be used for this procedure such as interference screws, taps, and drill bits and are provided non-sterile and must be properly cleaned and sterilized before each use.

### 10. Material Specifications

The OrthoSta tape is made of ultra-high molecular weight polyethylene (UHMWPE) coated in Healix silver ion antimicrobial. Interference Screws and buttons are made of titanium alloy (Ti 6Al-4V).

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