KNOTLESS/CRIMPLESS LATERAL STABILIZATION SYSTEM

Instructions for Use
Gebrauchsanweisung
Mode d'emploi
Modo de empleo
Istruzioni per l'uso

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ENGLISH

1. Device Description

The OrthoZip system consists of two sizes of Ultra-High Molecular Weight Polyethylene (UHMWPE) self-locking sutures coated in SilverPlus, a silver ion antimicrobial. The OrthoZip system is used in conjunction with titanium alloy bone anchors and surgical accessories including a screwdriver and drill bit.

2. Indications

The OrthoZip System is intended for veterinary use only in the management of ligament injuries. Applications include stabilization of cranial cruciate injuries and collateral ligament injuries or augmentation.

3. Contraindications

Contraindications include insufficient quality or quantity of bone, blood supply issues or previous infections that will prevent healing, foreign body sensitivity, infection, and or conditions that may prevent the patient from following post-op instructions. Do not use for surgeries other than those indicated.

4. Adverse Effects

Infection, foreign body reactions, soft tissue inflammation, aseptic loosening, and failure/breakage of the implants may occur.

5. Warnings

The OrthoZip system is intended to be used by a trained medical professional and is for veterinary use only. The UHMWPE self-locking suture is a single use device and must not be re-used or re-sterilized. Failure to use the proper drill bit diameter may cause failure or breakage of the implant system. Reference the product surgical technique prior to surgery to ensure proper sequencing of applying the self-locking suture to avoid loosening. The surgeon must be aware that if the implant needs to be removed it will require secondary surgery. Explanted devices should be safely disposed of and not re-used. Knowledge of surgical techniques and proper selection and placement of the implants are important considerations in the successful use of the OrthoZip System.

6. MRI Safety Information

The OrthoZip titanium alloy bone anchors have not been evaluated for safety and compatibility in the MR environment. Suture manufactured from UHMWPE is MR safe.

7. Precautions

Veterinarian medical professionals are advised to review the product surgical technique prior to surgery. Surgeons may apply their professional judgment in the selection of the appropriate suture size and bone anchor diameter. Ensure the indicated drill bit is used for the corresponding bone anchor. When placing the self-locking suture center the sheath on the femoral bone anchor and then pull the free ends distally to desired tension. If pulling free ends with instruments be sure to not overtighten. During tensioning be sure to check cranial drawer is eliminated in all range of motion. Once the self-locking suture is secured, trim the free ends leaving no more than 1-2mm remaining before closing.

8. Package and Labeling

Devices should only be accepted if the packaging and labeling

arrive intact. Please contact Movora customer service if the package has been opened or altered upon arrival. Do not use any device if the package or implant is damaged.

9. Sterilization

The OrthoZip UHMWPE self-locking suture is provided sterile using Ethylene Oxide. Check the labeling for more information and expiration date. This UHMWPE self-locking suture is single use only; do not re-sterilize. Certain instruments and accessories used for this procedure such as bone anchors, drill bits, and screwdrivers are provided non-sterile and must be properly cleaned and sterilized before each use.

10. Material Specifications

The OrthoZip self-locking suture is made of ultra-high molecular weight polyethylene (UHMWPE) coated in SilverPlus silver ion antimicrobial. The bone anchors are made of titanium alloy (Ti 6Al-4V).

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