



I-Loc® 3-4-5 CIS-TRANS Correlation Chart

#3 Nail CIS-TRANS Correlation Chart

CIS DIMENSION	BONE DIAMETER (MM)																
	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18
2		2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
3			2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
4				2	3	4	5	6	7	8	9	10	11	12	13	14	15
5					2	3	4	5	6	7	8	9	10	11	12	13	14
6						2	3	4	5	6	7	8	9	10	11	12	13
7							2	3	4	5	6	7	8	9	10	11	12
8								2	3	4	5	6	7	8	9	10	11
9									2	3	4	5	6	7	8	9	10
10										2	3	4	5	6	7	8	9
11											2	3	4	5	6	7	8
12												2	3	4	5	6	7
13													2	3	4	5	6
14														2	3	4	5
15															2	3	4
16																2	3

#4 Nail CIS-TRANS Correlation Chart

CIS DIMENSION	BONE DIAMETER (MM)															
	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17
2		2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
3			2	3	4	5	6	7	8	9	10	11	12	13	14	15
4				2	3	4	5	6	7	8	9	10	11	12	13	14
5					2	3	4	5	6	7	8	9	10	11	12	13
6						2	3	4	5	6	7	8	9	10	11	12
7							2	3	4	5	6	7	8	9	10	11
8								2	3	4	5	6	7	8	9	10
9									2	3	4	5	6	7	8	9
10										2	3	4	5	6	7	8
11											2	3	4	5	6	7
12												2	3	4	5	6
13													2	3	4	5
14														2	3	4
15															2	3



I-Loc® 3-4-5 CIS-TRANS Correlation Chart

#5 Nail CIS-TRANS Correlation Chart

CIS DIMENSION	BONE DIAMETER (MM)														
	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20
1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
2		2	3	4	5	6	7	8	9	10	11	12	13	14	15
3			2	3	4	5	6	7	8	9	10	11	12	13	14
4				2	3	4	5	6	7	8	9	10	11	12	13
5					2	3	4	5	6	7	8	9	10	11	12
6						2	3	4	5	6	7	8	9	10	11
7							2	3	4	5	6	7	8	9	10
8								2	3	4	5	6	7	8	9
9									2	3	4	5	6	7	8
10										2	3	4	5	6	7
11											2	3	4	5	6
12												2	3	4	5
13													2	3	4
14														2	3
15															2

Example

CIS DIMENSION	BONE DIAMETER (MM)										
	4	5	6	7	8	9	10	11	12	13	14
1	2	3	4	5	6	7	8	9	10	11	12
2		2	3	4	5	6	7	8	9	10	11
3			2	3	4	5	6	7	8	9	10
4				2	3	4	5	6	7	8	9
5					2	3	4	5	6	7	8

Resultant TRANS Dimension



9 Whippany Road, Suite B2-7 | Whippany, New Jersey 07981 | telephone 973.331.7800

biomedtrix.com

MKT-TD-0046-A