



# VETERINARY ORTHOPEDIC IMPLANTS

Over 25 Years  
of delivering the  
widest selection of  
veterinary orthopedic  
implants, power equipment  
and instruments  
with unparalleled  
customer service.

FROM CUT TO CLOSE™

Selection • Value • Education • Unparalleled Service

## SUTURE SELECTIONS

V.O.I. distributes wound closure products from a leading U.S. company that produces suture for both the veterinary and human market.

- 100% SATISFACTION GUARANTEE\*
- Exclusively distributed by VOI
- Proudly made in the USA
- Minimal tissue reaction
- Maximum tensile strength
- Low capillarity
- Smoothness-minimal tissue drag
- Plasticity and elasticity
- Ease of handling - low memory
- Knot security
- Cost effective

\*If you're not satisfied with the suture for any reason, we will provide a credit, replace the product, or issue your money back.

Did you know?  
All suture boxes and  
sachets, are labeled in the  
same manner. Let our  
Customer Account Specialists  
help you move to Demetech!

**CALL  
TODAY!**



### ETHICON TRADEMARK NAME:

PDS®  
VICRYL®  
MONOCRYL®  
PROLENE®

### VOI EQUIVALENT

DemeDIOX™  
DemeSORB™  
DemeCAPRONE™  
DemeLENE™



VOI is an exclusive distributor of Demetech® in the Veterinary Field

USA • CANADA • EUROPE • JAPAN

[www.vetimplants.com](http://www.vetimplants.com)

**904.436.6540**

email: [sales@vetimplants.com](mailto:sales@vetimplants.com)



# FAQs

## ANSWERS TO YOUR MOST ASKED SUTURE QUESTIONS

Did you know? All suture boxes and sachets, no matter the brand or distributor, are labeled the same way. There is a standard set for seamless use and transition.

### 2. What if the suture materials are different than what I am currently using?

DemeTECH® uses the same suture materials as the industry leading brand.

The following is an excerpt from DemeTECH®'s quality certification:

“Physical properties and functionality testing assured that the device conformed with the suture diameter, suture length, knot pull tensile strength, needle attachment strength, extractable color, and sterility to methods outlined in the USP 35.”

*“So far so good! We have been using the VOI lines of suture and have had zero problems. Handling and functionality seem comparable to more expensive brands. Thanks VOI!”*

Dr. Levien, Australia  
Using PDO, PGCL, PGA, Nylon

*“I really like the suture! It doesn't leak dye, the needles are sharp, and it handles great.”*

Dr. Vaughan, Texas - USA  
Using PDO

### 1. Is it easy to convert codes?

Definitely!

There are multiple ways we can convert codes for you:

1. Send a picture of the front of your suture box that displays the suture specifications.
2. Tell us the suture type, gauge size, needle style, needle length, and needle curvature.
3. Tell us your current vendor's product codes.

*Conversion charts available upon request.  
We will provide a personalized spreadsheet that displays your current vendor's code with our equivalent.*

### 3. What if I receive discounted pricing from my current vendor?

We want to know!  
Everyone loves getting a great price!

Our suture is priced cost effectively without jeopardizing quality standards. If you receive discounted pricing lower than our list price, we will price match with proof of invoice.